

BKC-6C Benchtop low speed centrifuge



Model: **BKC-6C**

Features:

Brushless DC motor, free maintenance, no powder pollution, micro computer control, digital display indicates the time, speed and RCF. Electric lid lock, over speed and imbalance protection. The centrifuge body is made of high quality steel, Safe and reliable.

Technical Parameter:

max speed:	6000r/min	max RCF:	3660×g
max volume:	6×50ml	power supply:	AC 220V 50HZ 2A
timer:	0~99min	noise:	≤45dBA
dimensions:	485×320×255	net weight:	23 KG
speed accuracy:	±10r/min		

Available fixed rotors:

NO	Rotors type	Max speed(r/min)	Volume(ml)	Max RCF(×g)
NO31601	Angle Rotor	4000	12×10/7/5ml	2150
NO31602		5000	12×15/7/5ml	3080
NO31603		4000	12×20ml	2220
NO31604		4000	18×10ml	2250
NO31605		4000	30×7/5ml	2250
NO31606		4000	6×50ml	2100
NO31607		4000	24×10ml	2200
NO30615		6000	6×10/7/5ml	3450
NO30638		6000	6×15/7/5ml	3660
NO32665	Swing Rotor	4000	6×10/7/5ml (vacuum tube)	2390

